
## HemoCue Hb 201+

Operating Manual

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**HemoCue Hb 201+ *Operating Manual***  
**HemoCue Hb 201+ *Bedienungsanleitung***  
**HemoCue Hb 201+ *Manuel d'utilisation***  
**HemoCue Hb 201+ *Gebruiksaanwijzing***

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**Components**  
**Bestandteile**  
**Composants**  
**Onderdelen**

**GB**

1. The HemoCue Hb 201+ Analyzer.\*
2. Mains adapter.\*\*
3. 4 type AA or R6 batteries, 1.5 V.\*\*\*
4. A vial of HemoCue Hb 201 Microcuvettes.\*\*\*  
Individually packaged HemoCue Hb 201 Microcuvettes.\*\*\*
5. HemoCue Hb 201+ Operating Manual, HemoCue Hb 201+ Quick Reference Guide and Instruction CD.
6. HemoCue Cleaner.

The HemoCue Hb 201+ and its mains adapter are delivered in a carton. Open the carton on a stable surface and lift out the instrument and accessories. Consult local environmental authorities for proper disposal of batteries.

\* Do not open the cover of the analyser.

\*\* ⚠ Only use adapters listed under specifications.

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
**GB**

- 1a. If mains power is available, connect the adapter to the socket on the back.
- 1b. If no mains power is available, insert the 4 type AA or R6 batteries, 1.5 V, into the battery compartment.
- 1c. If a battery symbol appears on the display the batteries are running low. The analyser will continue to give accurate results but the batteries should be replaced as soon as possible.
2. The analyser can be connected to a PC (for further information, see page 9 or separate manual for PC connection) or directly to a printer (see page 26).
3. Pull the cuvette holder out to its loading position.
4. Press and hold the left button until the display is activated (all symbols appear on the display).
5. The display shows the version number of the programme, after which it will show "Σ", "Hb" and "🔊" (audio signal). During this time the analyser will automatically verify the performance of the optronic unit.
6. After 10 seconds the display will show three flashing dashes and the HemoCue symbol. This indicates that the HemoCue Hb 201+ is ready for use.

To turn off the analyser press and hold the left button until the display reads OFF and becomes blank.

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
**GB**

The time and date function within the set-up menu must be activated before data can be transferred to a PC.

1. Press both buttons at the same time.
2. The display now shows a flashing QC-symbol.
- 3ab. Use the right button to scroll until the display shows an audio symbol in the upper right corner. The signal can be activated or deactivated by pressing the left button.
4. Continue to scroll using the right button until the display shows characters for time, date and year. The hour figure will be flashing.
5. Use the right button to change the hours, minutes, day, month and year. To change a flashing figure hold the left button down for quick advancement or advance slowly by quickly pressing and releasing the button.
6. When all the settings are completed, hold the right button down for approximately 5 seconds. The analyser will then automatically return to the measuring position.

The analyser is supplied with default time and date format. This format can only be changed by using the PC programme (for further information, see separate manual for PC connection).

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
**Meting** Capillair bloed

**GB**

1. After start-up, the cuvette holder should be in its loading position. The display will show three flashing dashes and the HemoCue symbol.
2. Make sure the patient's hand is warm and relaxed. Use only the middle or ring finger for sampling. Avoid fingers with rings on.
3. Clean with disinfectant and allow to dry.
4. Using your thumb, lightly press the finger from the top of the knuckle towards the tip. This stimulates the blood flow towards the sampling point.

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


**GB**

5. For best blood flow and least pain, sample at the side of the fingertip, not in the centre.
6. Whilst lightly pressing towards the fingertip, prick the finger using a lancet.
7. Wipe away the first 2 or 3 drops of blood.
8. Re-apply light pressure towards the fingertip until another drop of blood appears

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


GB

9. When the blood drop is large enough, fill the microcuvette in one continuous process. Do NOT refill!
10. Wipe off excess blood on the outside of the microcuvette tip. Make sure that no blood is drawn out of the microcuvette during this procedure.
11. Look for air bubbles in the filled microcuvette. If present, take a new sample. Small bubbles around the edge can be ignored.

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
**GB**

- Place the filled microcuvette in the cuvette holder. This should be performed within ten minutes after filling the microcuvette!
- Push the cuvette holder to its measuring position.
- During the measurement “ $\frac{\infty}{\infty}$ ” will be shown on the display.
- After 15–60 seconds the haemoglobin value of the sample is displayed. The result will remain on the display as long as the cuvette holder is in the measuring position. When operating on battery power the analyser will automatically turn off after approximately 5 minutes.
- Although the reagents are present in the microcuvette in extremely low quantities, consult local environmental authorities for proper disposal. Always handle blood specimens with care, as they might be infectious.

All results are stored in the analyser and can be viewed in the display using the scroll function (see page 22) or by using the PC programme (for further information, see separate manual for PC connection).

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**Measuring** Venous or arterial blood

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


4

**GB**

1. Mix the blood well before performing the measurement.
- 2ab. Place a drop of blood onto a hydrophobic surface, e.g., a plastic film, using a pipette or a DIFF-Safe™.
3. Fill the microcuvette in one continuous process. Do NOT refill! Wipe off excess blood on the outside of the microcuvette tip. Make sure that no blood is drawn out of the microcuvette during this procedure.
4. Look for air bubbles in the filled microcuvette. If present, take a new sample. Small bubbles around the edge can be ignored.

Perform the analysis as per steps 12–16 on page 17.

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
**Set-up** QC-test

**GB**

The HemoCue system can be controlled daily with whole blood, haemolysate or other controls. To obtain reliable results only use controls recommended by HemoCue.

1. Press both buttons at the same time.
2. The display shows a flashing QC-symbol.
3. Select QC-test by pressing the left button.
4. The analyser automatically returns to the measuring position and the QC-symbol appears in the display.
5. Fill the microcuvette with the control solution recommended by HemoCue and carry out the measurement in accordance with steps 12–16 on page 17.
6. Once the measurement has been made the analyser automatically returns to the measuring position and the QC-symbol disappears from the display. Patient sampling can begin.

If more QC-tests are to be performed, proceed in the same way following steps 1–6. To deactivate the QC-test, follow steps 1–2 and scroll using the right button until another set-up activity is shown on the display. Hold the right button down for approximately 5 seconds. The analyser automatically returns to its measuring position and the QC-symbol is no longer visible.

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
**Set-up** Geheugenfunctie – Scrollen

GB

The analyser automatically stores up to 600 results. When the memory is full the analyser will automatically overwrite the oldest result. To view results stored in the memory use the following procedure.

- 1ab. The cuvette holder may be in the loading position or in the measuring position. If in the measuring position the most recent measurement is displayed.
2. Press the left or right button to scroll backwards or forwards between the results. The stored values are now shown in the display. When the button is held down the user can see the order of the results (the most recent measurement is always number 1) and when the button is released the result is shown.
- 3ab. If no buttons are pressed whilst in the loading position the analyser automatically resets after 5 seconds and new measurements can be made. If in the measuring position pull out the cuvette holder and wait for the HemoCue symbol after which a new measurement can be made.

Note: The analyser will not respond to other commands while a measurement is being made. If the error code E35 is shown while scrolling through the results, the measurement has not been stored in the memory, see the trouble shooting guide on page 30.

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
**Set-up** Resultaten wissen

GB

1. Press both buttons at the same time.
2. The display now shows a flashing QC-symbol.
3. Scroll using the right button until the display shows a flashing waste bin in the left corner.
- 4ab. To delete the most recent result press and hold the left button. A warning triangle together with a "1" appears in the display for approximately 3 seconds. The most recent result has now been deleted and the display shows a "0" to confirm that this has been done. The most recent result can only be deleted immediately after the measurement has been made.
- 5ab. To delete all stored results hold both buttons down at the same time for 5 seconds. The display shows the waste bin, the warning triangle and the total number of results stored. Continue to hold both buttons down until the number changes to "000" and the warning triangle disappears. All stored results have now been deleted.

Note: Deleted results cannot be recovered. If results are to be saved they must be printed out or downloaded to a PC before starting the deleting process.

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
**Configurer** le mode impression

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**GB**

The analyser can be on or off when connected to the printer. Only the most recent result can be printed out directly after the measurement.

1. Connect the cable\* to the analyser and printer\*.
2. Perform the analysis by following the steps 12–16 on page 17.
3. When the result is shown on the display, the printer will automatically print the result, date and time.

\* Not included.

Note: Only use ASCII printers recommended by HemoCue.

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
Maintenance  
Wartung  
Maintenance  
Onderhoud

**GB**

1. Check that the analyser is turned off. The display should be blank.
2. Pull the cuvette holder out to its loading position. Use a pointed object to carefully depress the small catch positioned in the upper right corner of the cuvette holder.
3. Whilst keeping the catch depressed, carefully pull the cuvette holder in the direction in which the handle of the cuvette holder is pointing.
4. Clean the cuvette holder with alcohol or mild detergent.
5. Push the HemoCue Cleaner swab into the opening of the cuvette holder. Pull out and push in 5–10 times. If the swab is stained, repeat with a new swab. Stop when the swab comes out unstained.
6. Wait 15 minutes before re-using the analyser. Replace the cuvette holder.

The cover may be cleaned with alcohol or a mild soap solution.

One HemoCue Cleaner is included with the analyser. Contact your local HemoCue distributor to order more.

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# **GB** Trouble Shooting Guide

If you are unable to resolve the problem by following this Trouble Shooting Guide, please contact your local HemoCue distributor or HemoCue AB.

<b>Symptom</b>	<b>Explanation</b>	<b>Action</b>
The analyser shows an error code.	May be a temporary fault.	Turn off the analyser and turn it on again after 30 seconds. Take a new microcuvette and repeat the measurement. If the problem continues, see specific error code below.
E00	No stable endpoint is found within the time range. 1. The cuvette is faulty. 2. The circuit board is out of order.	1a. Check the expiry date for the microcuvettes. 1b. Take a new microcuvette and repeat the measurement. 2. The analyser needs service. Contact your distributor.
E01–E05	1. Fault or dirt in the optronic or electronic units.	1a. Turn off the analyser and clean the optronic unit using the HemoCue Cleaner. 1b. The analyser needs service. Contact your distributor.
E06	1. Unstable blank value. The analyser might be cold.	1. Turn off the analyser and allow it to reach room temperature. If the problem continues, the analyser needs service. Contact your distributor.
E07	1. The battery power is too low.	1a. The batteries need to be replaced. Turn off the analyser and replace the batteries, 4 type R6 or AA, 1.5 V. 1b. Use the mains power supply.
E08	The absorbance is too high. 1. Light blocking item in the cuvette holder.	1a. Check that the analyser and microcuvettes are being used according to the HemoCue Hb 201+ operating manual and instructions for use. 1b. The analyser needs service. Contact your distributor.
E09–E30	1. Fault or dirt in the optronic or electronic units.	1a. Turn off the analyser and clean the optronic unit using the HemoCue Cleaner. 1b. The analyser needs service. Contact your distributor.
E35	1. The current measurement has been incorrectly stored and will remain as error code E35 until the entire memory has been deleted.	1a. No action is normally required. 1b. If storage problems continue, the analyser needs service. Contact your distributor.

Symptom	Explanation	Action
HHH	1. Measured value exceeds 256 g/L (25.6 g/dL, 15.9 mmol/L).	
No characters on the display.	<ol style="list-style-type: none"> <li>1. The analyser is not receiving power.</li> <li>2. If on battery power, the batteries need to be replaced.</li> <li>3. The display is out of order.</li> </ol>	<ol style="list-style-type: none"> <li>1a. Check that the mains adapter is connected to the mains power supply.</li> <li>1b. Check that the mains adapter is securely connected to the analyser.</li> <li>1c. Check that the cable is not damaged.</li> <li>2. Turn off the analyser and replace the batteries, 4 type R6 or AA, 1.5 V.</li> <li>3. The analyser needs service. Contact your distributor.</li> </ol>
The display gives erroneous characters.	<ol style="list-style-type: none"> <li>1. The display is out of order.</li> <li>2. The microprocessor is out of order.</li> </ol>	<ol style="list-style-type: none"> <li>1. The analyser needs service. Contact your distributor.</li> <li>2. The analyser needs service. Contact your distributor.</li> </ol>
The display shows "▣".	<ol style="list-style-type: none"> <li>1. The batteries need to be replaced.</li> <li>2. If on mains power, the mains adapter or the circuit board is out of order.</li> </ol>	<ol style="list-style-type: none"> <li>1. Turn off the analyser and replace the batteries, 4 type R6 or AA, 1.5 V.</li> <li>2a. Check that the mains adapter is properly connected and working.</li> <li>2b. The analyser needs service. Contact your distributor.</li> </ol>
The display does not switch from "Σ" and "Hb" to three flashing dashes and "☺" (ready for measuring).	<ol style="list-style-type: none"> <li>1. The magnet in the cuvette holder may be missing.</li> <li>2. The magnetic sensor is out of order.</li> </ol>	<ol style="list-style-type: none"> <li>1. The analyser needs service. Contact your distributor.</li> <li>2. The analyser needs service. Contact your distributor.</li> </ol>

Symptom	Explanation	Action
<p>Measurements on control materials are out of range – either too HIGH or too LOW.</p>	<ol style="list-style-type: none"> <li>1. The microcuvettes are beyond their expiry date, damaged or have been improperly stored.</li> <li>2. The optical eye of the microcuvette is contaminated.</li> <li>3. The control has not been mixed properly and/or is not at room temperature.</li> <li>4. Air bubbles in the microcuvette.</li> <li>5. The optronic unit is dirty.</li> <li>6. The control is not suitable for use with the HemoCue Hb 201+ system.</li> <li>7. The calibration of the analyser has been changed.</li> <li>8. The controls are beyond their expiry dates or have been improperly stored.</li> </ol>	<ol style="list-style-type: none"> <li>1. Check the expiry date and the storage conditions of the microcuvettes.</li> <li>2. Remeasure the sample with a new microcuvette.</li> <li>3. Make sure that the control is mixed properly and at room temperature.</li> <li>4. Check the microcuvette for air bubbles. Remeasure the sample with a new microcuvette.</li> <li>5. Clean the optronic unit using the HemoCue Cleaner.</li> <li>6. Only use controls intended for the HemoCue Hb 201+ system recommended by HemoCue.</li> <li>7. The analyser needs service. Contact your distributor.</li> <li>8. Check the expiry date and the storage conditions of the control. Take a new microcuvette and repeat the measurement.</li> </ol>
<p>Measurements on patient samples are higher or lower than anticipated.</p>	<ol style="list-style-type: none"> <li>1. The microcuvettes are beyond their expiry date, damaged or have been improperly stored.</li> <li>2. The optical eye of the microcuvette is contaminated.</li> <li>3. Air bubbles in the microcuvette.</li> <li>4. The optronic unit is dirty.</li> <li>5. The calibration of the analyser has been changed.</li> </ol>	<ol style="list-style-type: none"> <li>1. Check the expiry date and the storage conditions of the microcuvettes.</li> <li>2. Remeasure the sample with a new microcuvette.</li> <li>3. Check the microcuvette for air bubbles. Remeasure the sample with a new microcuvette.</li> <li>4. Clean the optronic unit using the HemoCue Cleaner.</li> <li>5. The analyser needs service. Contact your distributor.</li> </ol>

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

# Specifications

## General

The HemoCue Hb 201 is a system used for the determination of the total amount of haemoglobin in whole blood. The system consists of a specially designed analyser with specially designed microcuvettes containing dried reagents. The microcuvette serves as pipette, reaction vessel and as a measuring microcuvette. No dilution is required. The haemoglobin measurement takes place in the analyser, which follows the progress of the reaction until the end point has been reached. The system is factory calibrated against the haemiglobincyanide (HiCN) method, the international reference method for the determination of the haemoglobin concentration in blood, and needs no further calibration.

## Intended purpose

Quantitative determination of haemoglobin in capillary, venous or arterial blood. The HemoCue Hb 201 system is used for quantitative determination of haemoglobin in whole blood using a specially designed analyser, the HemoCue Hb 201+ Analyzer, and specially designed microcuvettes, the HemoCue Hb 201 Microcuvettes. HemoCue Hb 201 Microcuvettes are for *In Vitro* Diagnostic use only. The HemoCue Hb 201+ Analyzer is only to be used with HemoCue Hb 201 Microcuvettes.

## IVD Medical Device Directive

The HemoCue Hb 201+ complies with the IVD Medical Device Directive 98/79/EC and carries the CE mark.

## Theory

Sodium deoxycholate haemolyses the erythrocytes and haemoglobin is released. Sodium nitrite converts haemoglobin to methaemoglobin which, together with sodium azide, gives azidemethaemoglobin. The absorbance is measured at two wavelengths (570 nm and 880 nm) in order to compensate for turbidity in the sample.

## Storage and handling

Use the HemoCue Hb 201 Microcuvettes prior to their expiry date. The expiry date is printed on each package.

### Storage for microcuvettes kept in a vial

The microcuvettes are to be stored at room temperature (15–30 °C, 59–86 °F). Do not refrigerate.

Once the seal is broken the microcuvettes are stable for three months. Always keep the container properly closed.

### Storage for individually packaged cuvettes

The microcuvettes are to be stored at room temperature (15–30 °C, 59–86 °F). Do not refrigerate.

### HemoCue Hb 201+ Analyzer

The analyser can be stored at 0–50 °C (32–122 °F). Operating temperature is 15–30 °C (59–86 °F). Allow the analyser to reach ambient temperature before use. The analyser should not be operated at high (i.e. > 90 % non-condensing) humidity.

## Specimen collection and preparation

Capillary, venous or arterial blood may be used. Appropriate anticoagulants (e.g. EDTA or heparin) may be used, preferably in solid form to avoid dilutional effects. Mix all specimen tubes thoroughly on a mechanical mixer for at least 2 minutes or invert the tube 8–10 times by hand. Hemoglobin remains unchanged for days, provided that the blood does not become infected. If the specimen has been stored in a refrigerator, it will be viscous and the blood should be allowed to warm up to room temperature before mixing<sup>2</sup>.

In case of peripheral circulatory failure, hemoglobin measurements from capillary sampling may be misleading. Confirmation of an unacceptable or unexpected laboratory result is considered as a common, well established routine in the laboratory. Extra cellu-

lar fluid present in capillary samples due to individual physiology or the affects of unacceptable capillary sampling technique can cause false hemoglobin results.

One of several pre analytical factors associated with technique revolves around the size and action mode of the lancet and how it is used. This and other factors associated with capillary sampling technique can have an impact on the puncture, blood flow and ultimately on the results obtained. Additional factors unrelated to sampling technique, such as an individual's capillary circulation status, must also be recognized as potential causes of false results. Therefore when testing for hemoglobin using capillary blood as the sample, confirmation of an unacceptable or unexpected result is needed to rule out potential pre analytical factors as the cause. This is particularly helpful when a cut off value is used, i.e. when screening for nutritional status based on hemoglobin. Since the individual performing the test may not always be aware of all the potential pre analytical factors or even those that are directly or indirectly associated with sampling techniques, a second finger stick should be performed when wanting to rule out pre analytical factors as a cause of an unacceptable or unexpected result. Performance of the second finger stick by a different individual may aid in lessening the affect of the factors that could have been technique related.

**Always handle blood specimens with care, as they might be infectious. Consult local environmental authorities for proper disposal. Always wear protective gloves when handling blood specimens. The microcuvette is for single use only.**

#### **Materials required**

- HemoCue Hb 201<sup>+</sup> Analyzer
- HemoCue Hb 201 Microcuvettes
- Lancet (for capillary samples)

- Pipette or other transfer device (for venous, arterial or control material samples)
  - Lint-free wipe (non-fraying)
- Hydrophobic plastic or glass slide (for venous, arterial or control material samples)

#### **Quality Control**

The HemoCue Hb 201<sup>+</sup> Analyzer has an internal quality control, the "self test". Every time the analyser is turned on, it will automatically verify the performance of the optronic unit of the analyser. This test is performed every second hour if the analyser remains switched on.

If quality control checks are required for regulatory reasons they should be performed using liquid controls recommended by HemoCue.

#### **Expected values<sup>2</sup>**

Adult Males 130–170 g/L (13.0–17.0 g/dL, 8.1–10.5 mmol/L)

Adult Females 120–150 g/L (12.0–15.0 g/dL, 7.4–9.3 mmol/L)

Infants, after neonatal period 110–140 g/L (11.0–14.0 g/dL, 6.8–8.7 mmol/L)

Children, two years to teenage: gradual increase to adult normals. Due to a wide range of conditions (dietary, geographical, etc) which affect normal values, it is recommended that each laboratory establish its own normal range.

#### **Measuring range**

0–256 g/L (0–25.6 g/dL, 0–15.9 mmol/L).

Results above 256 g/L (25.6 g/dL, 15.9 mmol/L) will be displayed as HHH.

#### **Limitations**

HemoCue Hb 201 Microcuvettes are for *In Vitro* Diagnostic use only. The HemoCue Hb 201<sup>+</sup> Analyzer is only to be used together

with HemoCue Hb 201 Microcuvettes. For further limitations of the procedure, see the HemoCue Hb 201 Microcuvettes instructions for use.

HemoCue Lancets  
HemoCue Hb 201<sup>+</sup> Control Cuvette

### **Technical Specifications**

Dimensions: 85×160×43 mm

Weight: 350 g (batteries included)

4 batteries type AA or R6, 1.5 V

Mains adapter: CE marked

Only use adapters as listed on page 50.

Pollution degree: 2

Overvoltage category: II

The instrument is tested according to IEC 61010-1, First edition, 1990 with Amendment 1, 1992 and Amendment 2, 1995 (EN 61010-1: 1993 with A2: 1995), EN 60601-1-2 and complies with the IVD Medical Device Directive 98/79/EC.

### **Warranty**

The analyser carries a 24-month warranty from the day of receipt. After the warranty period service/repair is carried out at fixed prices.

### **Service and Disposal**

The analyser should be cleaned as recommended under Maintenance prior to service or disposal. Consult local environmental authorities for proper disposal.

### **Spare parts – Accessories**

The following spare parts and accessories are available:

Mains adapter

Battery Lid

PC programme

Cables for PC and printer connection

Liquid controls

HemoCue Cleaners

Hard Carrying Case

**Recommended separation distance between Portable and mobile RF communications equipment and HemoCue Hb 201<sup>+</sup> Analyzer**

The HemoCue systems are intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of HemoCue systems can help prevent electromagnetic interference by maintaining a minimum distance between portable and RF communications equipment (transmitters) and HemoCue systems as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz $d=1.2\sqrt{P}$	80 MHz to 800 MHz $d=1.2\sqrt{P}$	800 MHz to 2.5 GHz $d=2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at maximum output power not listed above, the recommended separation distances (d) in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where (P) is the maximum output power rating of the transmitter in watts (W) according to transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

### Guidance and manufacturer's declaration – Electromagnetic immunity

The HemoCue systems are intended for use in the electromagnetic environment specified below. The customer or user of the HemoCue systems should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/ burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	2 kV for power supply lines ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode	±1 kV differential mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply lines IEC 61000-4-11	<5 %U (>95 % dip in U) for a 0.5 cycle 40 %U (60 % dip in U) for 5 cycles 70 %U (30 % dip in U) for 25 cycles <5 %U (>95 % dip in U) for 5 seconds For explanation of U see NOTE 1	<5 %U (>95 % dip in U) for a 0.5 cycle 40 %U (60 % dip in U) for 5 cycles 70 %U (30 % dip in U) for 25 cycles <5 %U (>95 % dip in U) for 5 seconds	Mains power quality should be that of a typical commercial or hospital environment.  If the user of the HemoCue systems requires continued operation during power mains interruptions, it is recommended that the HemoCue systems be powered from an uninterruptible power supply or a battery.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
<p>Conducted RF IEC 61000-4-6</p> <p>Radiated RF IEC 61000-4-3</p>	<p>3 Vrms 150 kHz to 80 MHz</p> <p>3V/m 80 MHz to 2.5 GHz</p> <p>See NOTE 2 and NOTE 3</p>	<p>3 Vrms</p> <p>3 V/m</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the HemoCue systems, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter</p> <p>Recommended separation distance</p> $d=1.2\sqrt{P}$ $d=1.2\sqrt{P} \text{ 80 MHz to 800 MHz}$ $d=2.3\sqrt{P} \text{ 800 MHz to 2.5 GHz}$ <p>Where (P) is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and (d) is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey (a), should be less than the compliance level in each frequency range (b).</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol.</p> 

NOTE 1 U is the a.c. mains voltage prior to application of the test level.

NOTE 2 At 80MHz and 800MHz, the separation distance for the higher frequency range applies.

NOTE 3 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the HemoCue systems are used exceeds the applicable RF compliance level above, the HemoCue systems should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the systems.
- b) Over the frequency range 150 KHz to 80 Mhz, field strength should be less than 3 V/m.

**Technical specifications (EMC-RF)**

**Use only cables with the following specification:**

USB shielded maximum 2 m

Serial shielded maximum 1.5 m

<b>Guidance and manufacturer's declaration – electromagnetic emissions</b>		
The HemoCue systems are intended for use in the electromagnetic environment specified below. The customer or the user of the HemoCue systems should assure that it is used in such an environment.		
Emission test	Compliance	Electromagnetic environment – guidance
RF emissions	Group 1	The HemoCue systems uses RF energy only for its internal function. Therefore, its RF emissions are very low and not likely to cause any interference in nearby electronic equipment
RF emissions	Class B	The HemoCue systems are suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

### Symbols Used

Attention, see instructions for use

CE mark

Class II equipment

Only valid within the European Community. Indicates separate collection for waste of electrical and electronic equipment.

Temperature limitation

Efficiency Level

### References

1. Makarem, A. In clinical Chemistry: Principles and Technics, 2nd ed., Henry, R.J., Cannon, D.C., and Winkelman, J.W., Harper and Row, Hagerstown, M.D., 1974, pp. 1125–1147.
2. Dacie and Lewis, Practical Haematology, Ninth Edition, 2001.
3. Reference and Selected Procedures for the Quantitative Determination of Hemoglobin in Blood; Approved Standard NCCLS Document H15-A.
4. Evaluation of precision performance of clinical chemistry devices; approved Guideline NCCLS Document EP5-A.

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HemoCue AB, Box 1204, SE-262 23 Ängelholm, Sweden. Phone +46 77 570 02 10. Fax +46 77 570 02 12.