

Features

Autonomously functioning cart with integrated oil-free compressor, silent suction machine and spray water tanks

High quality stainless steel housing

1 Three-way syringe (air/water/spray) with sterilisable tips

Integrated spray water system for syringe and instruments (check valve prevents the back flow of water)

3 Instrument hoses with light (6-hole connector)

Instrument hoses detachable

Spray adjustment regulators for each instrument

1 Saliva/Surgical suction hose including oral-dry-cup

1 High-performance suction hose 300 l/min

Continuously adjustable

Unique BPR Swiss technology

2 Spray water tanks: 2 × 0.75 litres

2 Suction fluid containers: 2 × 2 litres (automatic overflow protection)

Foot control with spray on/off function and chip

Integrated compressor: Capacity 75 l/min at 4 bar;

Sound pressure level: 48 dB; Compressed air tank: 4 litres

Technical Data

Air consumption instruments: 35–55 l/min

Operating pressure handpieces: 2.5–3.5 bar

Each instrument is individually adjustable

Total weight: 55 kg

Dimensions: 535mm × 580mm × 960mm

Warranty: 2 years

Options: Electric micromotor with light including display and 30 programs, electric scaler, peristaltic pump, turbine, air motor, curing light, tray etc.

Combi-Cart Clinic All-in-one






The COMBI-CART CLINIC combines 1 three-way-syringe, 3 treatment instruments and 2 suction devices including compressor and suction machine on one single cart!

Aside from treating your patients in your dental clinic, it offers you the flexibility to treat them in hospitals, retirement homes or other institutions. You can easily move your COMBI-CART CLINIC from one room to the next.

No installation needed (plug-and-play)! To activate and use your COMBI-CART CLINIC simply plug in the power cord. The silent COMBI-CART CLINIC, designed and manufactured in Switzerland, in addition to flexible application, is also marked by high performance, little need to apply maintenance, and easy use.

The high performance suction (300 l/min!) with unique BPR Swiss technology is complemented by instruments thoroughly efficient even with continuous operation. The great comfort of the COMBI-CART CLINIC is rounded off by the convenience of a very silent compressor (48 dB!) and useful connections for external devices.



-  Continuous Operation
-  Full Performance
-  Plug-and-Play
-  All-in-One
-  Mobile/Rollable