

Features

Autonomously functioning cart with integrated silent oil-free compressor and spray water tanks

High quality stainless steel housing

1 Three-way syringe (air/water/spray) with sterilisable tips

Integrated spray water system for syringe and instruments (check valve prevents the back flow of water)

2 Spray water tanks: 2 × 0.75 litres

4 Instrument hoses with light (6-hole connector)

Instrument hoses detachable

Spray adjustment regulators for each instrument

Foot control with spray on/off function and chip

Integrated compressor: Capacity 75 l/min at 4 bar;
Sound pressure level: 48 dB; Compressed air tank: 4 litres

Technical Data

Air consumption instruments: 35–55 l/min

Operating pressure handpieces: 2.5–3.5 bar

Each instrument is individually adjustable

Total weight: 53 kg

Dimensions: 535mm × 580mm × 960mm

Warranty: 2 years

Options: Electric micromotor with light including display and 30 programs, electric scaler, peristaltic pump, turbine, air motor, curing light, tray etc.

P3-Cart Pure freedom



The elegant P3-CART treatment cart requires no installation and is totally mobile. To use the P3-CART simply insert the power plug!

The instruments can be freely selected and are fully efficient even with continuous operation.

The high performance of the P3-CART is complemented by uncompromising comfort. The operation of the cart is silent (48 dB!) and oil-free.

The comfort is further enhanced by the great design and attention to detail. Examples of this are the external removable instrument hoses and useful connections for external devices. Further the P3-CART is equipped with two spray water tanks. If a spray water tank becomes empty, you can change to the second tank with a switch, continuing your treatment without interruption.



-  Continuous Operation
-  Full Performance
-  Plug-and-Play
-  Mobile/Rollable