



43 CLINICAL EVALUATORS
1304 TOTAL USES
92% CLINICAL RATING

Key features: Streamlined universal composite system • 3 shades match all 16 VITA classical shades • Designed for use without a blocker

Description

SimpliShade is a nano-hybrid composite resin that:

- Matches all 16 VITA classic shades with just 3 shades: Light, Medium, and Dark
- Is designed for use in all cavity classes without the need for a blocker
- Incorporates the Adaptive Response Technology (ART) filler network, also found in **Harmonize** which enables dynamic viscosity for superior handling and optical characteristics that mimic human enamel.

Indications

- Direct placement in anterior and posterior teeth
- Core build-ups
- Repair of porcelain and composite restorations

Unique Attributes

- Reduces composite inventory
- Makes shade selection simpler and faster
- Excellent esthetics
- Also available in a bleach shade and a universal opaque shade



BEFORE



AFTER

Clinician: Dr. Sabiha S. Bunek
Pre- and Post-treatment of tooth #8.



Clinical Tips


- If you are going to use this material consistently, arrange the VITA shade guide light to dark to correspond to SimpliShade light to dark equivalent shades.
- Polish with diamond paste and you will get a glassy smooth finish that rivals your most beautiful glass-ceramic restorations.
- If you have to choose between light and medium, go for the light because it tends to err on the side of darker.

“SUPERIOR HANDLING PROPERTIES.”


Evaluators' Comments

- “An outstanding composite that should be in every dentist's armamentarium of go-to restorative materials.”
- “Good product that is very useful in most clinical cases with minimal inventory or material.”
- “Good addition for difficult matching situations.”
- “Excellent handling characteristics with a superior chameleon effect, even without blockers.”
- “Too translucent for certain instances.”
- “A good product for anteriors; however, I prefer a bit more firmness in a posterior composite.”

Evaluation Summary:



Compared to Competitive Products:



Consultants who would:

- 86%** Recommend to a colleague
- 5%** Yes, instead of current product
- 60%** Yes, in addition to current product
- 16%** I might want to order for certain cases