



OBSERVA SERIES

PrimaLED

The first and only examination lamp with the possibility to choose an ambient light in addition to traditional concentrated light. The Ergo-Spring balancing system makes PRIMALED very easy to handle and stable. PRIMALED is ideal for any type of installation, from the outpatient department to the intensive care unit.

Fixing Systems





PRIMAFLEX

Gooseneck Arm

PRIMALED_P11_BRACCIO FLESSIBILE

PRIMA

Joint Arm

PRIMALED_P11_SNODI

Performances

Primaled

Light intensity at 0,5 m distance (Ec)	110.000 Lux
Out reflector diameter	19,5 cm
Color temperature: double selection	4000/4500 °K
Color rendering index (CRI)	95 Ra
d ₁₀ Light field diameter @ 50 cm (Ec)	150 mm
Total radiated energy Ee where the illuminance reaches max level	294 W/m ²
Ratio between radiated energy Ee and illuminance EcAverage led life	3,09 mW/m ² ·lx
Electrical absorption	22 VA
Control of the illuminance	25 - 100%