

Charisma® ABC
Economic simplicity.



CHARISMA®

Easy as ABC.

Charisma ABC is a new hybrid composite that allows natural restorations with minimal effort. Its shade selection, creamy consistency and easy polishing make Charisma ABC your go-to composite for every restoration.

Enjoy the following benefits with your simple yet trusted restorative solution from Heraeus Kulzer:

- Easy application
- A brand you can trust
- Savings you can count on

Giving a hand to oral health.

Charisma® ABC – Economic simplicity.

Charisma ABC is a light-curing, radiopaque hybrid composite for everyday restorations:

- **Easy application**
- **A brand you can trust**
- **Savings you can count on**

The formula behind Charisma ABC ensures:

- Easy handling for everyday restorations.
- Great polishability that saves time while providing a natural look.
- Simple shade range for easy shade selection and single-shade layering techniques.

Our expertise and long-term experience with composite technology ensure the quality and reliability you have come to expect from Heraeus Kulzer.

With a broad range of composite materials, we are committed to meeting all your needs – from simple fillings with Charisma ABC to the most complex and challenging restorations with our highly esthetic composites Venus Diamond® and Venus® Pearl. Whatever the situation, whatever your standard, we always aim to recognize and fulfill your unique needs.



Charisma ABC at a glance:	
Filler content	59% by Vol.
Filler size	0.005 – 10 µm
Filler system	Barium Aluminium Fluoride glass, highly dispersive silicon dioxide
Flexural strength [MPa] LED, 20s	119
Elasticity modulus [GPa] LED, 20s	8.0
Shrinkage [%-vol] acc. to Archimedes	2.8
Reflection [%] after toothbrush abrasion, 100.000 cycles	5.2
Roughness [µm] after toothbrush abrasion, 100.000 cycles	0.6
Radiopacity [%-Al]	196

Source: Heraeus Kulzer R&D, Internal measurements. Documentation available.



Before

After

Source: Dr. Janine Schweppe, Heraeus Kulzer GmbH, Germany

Charisma® ABC- Product range.

Shades	Charisma® ABC Syringe Refill 1 x 4g	Charisma® ABC PLT Refill 20 x 0.2g
A1	66064772	66064780
A2	66064773	66064791
A3	66064774	66064792
A3.5	66064775	66064793
A4	66064776	66064794
B2	66064777	66064795
C2	66064778	66064796
OPAQUE	66064779	66064797

