

A.S.A.P.[®] Indirect+

All Surface Access Polishers



Chairside Zirconia Polishing Made Easy


The **ASAP[®] Indirect+** diamond polishing system is all you need to refine and create an outstanding shine on all your indirect ceramic restorations.

- Designed with a high concentration of diamond particles coupled with a durable silicone binder for use with all substrates.
- Multi-shaped adjusters allow for refinement and the ability to remove scratches often caused by burs.
- A fast, super high gloss polish is completed with the diamond-impregnated spiral wheel polishers, Pre-polisher (blue) and Final Polisher (pink).

CLINICAL TECHNIQUE

Chairside Adjustment and Polishing with the A.S.A.P. Indirect+ System

Dentistry and photography courtesy of Robert Lowe, DDS



1 To make minor adjustments or to remove scratches caused by diamond bur adjustments, the A.S.A.P. Indirect+ Coarse Adjuster; available in a disc, point or cylinder shape, can be used on any zirconia, lithium disilicate or ceramic prosthesis.

3 The A.S.A.P. Indirect+ Pre-polisher is used to smooth the prosthetic surface and establish an initial gloss.

4 The A.S.A.P. Indirect+ Final Polisher creates a stunning high luster in seconds on any zirconia, lithium disilicate or ceramic surface.

A.S.A.P. INDIRECT+ INTRAORAL/RA STARTER KIT (130918)

Contains: 3 Coarse Adjusters (1 x Point, 1 x Disc, 1 x Cylinder), 1 x Pre-polisher, 1 x Final Polisher

A.S.A.P. INDIRECT+ EXTRAORAL/HP STARTER KIT (130940)

Contains: 2 Coarse Adjusters (1 x Disc, 1 x Cylinder), 2 x Pre-polishers, 2 x Final Polishers

A.S.A.P. INDIRECT+ INTRAORAL/RA POLISHER REFILL 2-PACK (130919)

Contains: 1 x Pre-polisher, 1 x Final Polisher

A.S.A.P. INDIRECT+ EXTRAORAL/HP POLISHER REFILL 2-PACK (130941)

Contains: 1 x Pre-polisher, 1 x Final Polisher

A.S.A.P. INDIRECT+ INTRAORAL/RA REFILL 3-PACKS

Pre-polisher (130920)
Final Polisher (130921)
Point Adjuster (130922)
Cylinder Adjuster (130923)
Disc Adjuster (130924)

A.S.A.P. INDIRECT+ EXTRAORAL/HP REFILL 3-PACKS

Pre-polisher (130942)
Final Polisher (130943)
Cylinder Adjuster (130944)
Disc Adjuster (130945)



To view the latest news, education and more, visit:

www.clinicianschoice.com

Clinician's Choice[®]
Designed for better dentistry.