




OPTIMA Apex Locator



Accurate

- Accuracy rate to 97.71%
- Multi-frequency apical positioning technology is capable of accurate measurement even in the presence of blood or residual pulp in the root canal
- Automatic calibration: ensures accuracy



Simple

- Clear bright LCD display
- Clear sound with intermittent beeping at different frequency depending on the position of the apex
- Sound level can be adjusted, and the position can be distinguished by sound, making measurement easier

Apical Stop Setting

- The apical stop alarm function position can be adjusted freely according to clinical needs
- Flexible setting, timely reminder to measure distance
- Reposition the display location of 0.5mm from the apical foramen

Contents included:

- 2 x file clip
- 2 x lip hook
- 1 x touch probe
- 1 x calibration tester
- All accessories are autoclavable

www.bainternational.com